

More beauty, plenty of brawn

BruxZir® Anterior... the new benchmark  
for maximum strength & esthetics



**BruxZir®**  
**ANTERIOR**  
— SOLID ZIRCONIA —



For more information,  
visit [www.bruxzir.com](http://www.bruxzir.com)



# BruxZir® Anterior, significantly more translucent for natural vitality

**BruxZir® Anterior** is the latest advancement in the BruxZir Solid Zirconia product line. This highly esthetic restorative material is designed specifically to satisfy the esthetic and functional requirements of the anterior region of the mouth. Exhibiting an average flexural strength of 650 MPa with translucency and color similar to natural dentition, BruxZir Anterior is an ideal, esthetic solution for your anterior cases.

- ✓ Anterior and premolar crowns
- ✓ Anterior and premolar bridges with one pontic
- ✓ More conservative preparation than lithium disilicate restorations; 0.8 mm required, 1.25 mm is ideal



## BruxZir® Anterior Preparation Guidelines



*These Captex™ Crowns had been in place for 15 years and the patient exhibited a possible metal allergy. Though the crowns had served her well, she was not pleased with the visible margins and was happy to hear that there are now several more esthetic options available.*

*The underlying preparations were tetracycline stained and would show through had IPS e.max® HT (High Translucency) ingots been used. The case was restored with BruxZir Anterior crowns, which combine an esthetic level of translucency with the ability to block out discolorations.*



*The crowns on tooth #6, #7 and #8 are BruxZir Anterior and the crowns on tooth #9, #10 and #11 are IPS e.max. Due to the increased translucency of BruxZir Anterior you can see that these solid zirconia crowns come much closer to matching the proven esthetics of IPS e.max.*